


ATTEND THE

WORLD EXCLUSIVE PREMIER
OF OUR NEWEST DOSING MACHINE

Hall: 10

Stand: B46



Lightweight design does the trick

Due to its light construction and small weighing cell it is able to measure small doses. Hence, it is 6x more accurate than market standards.

A belt instead of a dosing screw

Screws or similar create pulsating dosages and thus an inconsistent dosing behaviour. The Ferlin One uses a belt instead of a screw, which results in consistency.

An all-inclusive interface

Our user friendly interface contributes to an inclusive labour market. It is intuitive and self explaining.

Specifications

3kg Loadcell

7.5 liters Capacity hopper Total weight including control 15 kg

Compressed air 6 bar

240v 50/60Hz Power supply

45W **Energy consumption**

0-25 kg/h Dosing capacity

